Belifu

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Model:SM9126

ACTIVATE YOUR 24 MONTH WARRANTY & GET
TENS ELECTRODE PADS

Register within 2 weeks after receiving your tens unit

PLEASE VISIT NOW

→www.belifu.co/warranty

CUSTOMER SERVICE

 support@belifu.co

 @Belifu

Included in this standard package

- Device x 1pc
- Large pad x 2pcs
- Medium pad x 6pcs
- Small pad x 2pcs
- Electric wire x 3pcs


When use 4 line cable, make sure to use at least 1pc white(positive)+1pc blue line(negative)at the same time.

- USB cable x1pc
- Instruction manual x 1pc
- Bagx1pc

Features

- 2 output with 24 modes
- 20 level intensity
- 20 minutes default timer, can be adjustable from 10-60 minutes
- Rechargeable lithium battery
- Battery indicator to know when need to recharge battery
- 110-240v

Illustrations of the device and accessories

1. On/ Off switch
2. USB port
3. Battery indicator
4. Time
5. Mode selector
6. A Intensity adjustment (+)
7. Mode adjustment
8. A Intensity adjustment (-)
9. Time adjustment
10. B Intensity adjustment (+)
11. Mode adjustment
12. B Intensity adjustment (-)
13. Pause
14. A output connector
15. B output connector
16. Pad wires
17. Massage pad
18. USB cable

Modes

There are 24 modes to choose from. The “M” button cycles through the 6 main modes, all of these 6 main modes contain different built-in 4 modes.

1. Lateral palm kneading
2. Flat palm kneading
3. Crosscut fist kneading
4. Vertical fist kneading
5. Acupuncture
6. Moxibustion
7. Acupuncture kneading
8. Moxibustion kneading
9. Rubbing
10. Rubbing reduction
11. Tapping
12. Fibrillation massage
13. Lifting-pinching manipulation
14. Vibrating manipulation
15. High frequency manipulation
16. Fibrillation manipulation
17. Shiatsu-Deep
18. Thai Massage
19. Shiatsu
20. Strength Alternate Massage
21. Twinkling Cupping
22. Cupping
23. Fibrillation Cupping
24. Cup retaining

Operating instructions

1. Connect the pads to the wire by snapping them on, then remove the protective film.

Note: adding a few drops of water on the pads will increase their longevity.

2. Place pads where the massage sensations are desired. And plug into one of the ports on the device. For additional relief, use the second port and repeat steps 1&2 for two additional pads.

Make sure both pads on the skin and not overlapping.

Make sure your skin is free from any dirt, oil or lotion.

3. Turn the device on by pressing the on/ off switch on the top of the unit. The unit will be now activated and start in Model1.

4. Press the mode button to try different mode. Then press the center button to choose the waves within that selection. There are six icons on the screen, press “<” and “>” button to change the six icons.

Each icon contains different built-in 4 modes, press the “M” to choose the built-in 24 modes.

Please increase the intensity because it will return to lowest automatically when you change to a new mode.

5. Once you select a mode, gradually increase the A/B intensity by pressing the + button on the right/ left side; decrease it by pressing - button on the right/ left side of the tens.

- Adjust the tens using time by pressing the T button on the side of the tens. The timer can be increased by 10 minute intervals and the tens will automatically turn off after the set time has finished.
- Turn the device off by pressing the on/ off switch on the top of the unit If you need to turn off the device during massaging.

Precautions

- Do not apply the massage pads near the heart, on the head, above the neck, in the pubic region, over Scarred areas, on the throat or over the mouth.
- Avoid touching the pads when the unit is on.
- Pregnant women must avoid using this massager.
- Do not use this massager if you are using pacemaker or other life support equipment such as an artificial lung or respirator.
- Should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk or injury
- Do not use if you have following medical conditions:
Acute disease, trauma or surgical procedure in past six month
Cardiac dysfunction
Epilepsy

- High blood pressure
- High fever
- cancer
- Infectious disease
- Malignant tumors
- Recent surgery
- Skin problems
- Pregnancy

Maintenance and care

- When you are finished using the device, turn it off and take off the pads. Store the unit in a cool, dry place until next use.
- After usage, please cover the pads with the protective film. Applying a few drops water before and after usage will increase the longevity of the pads.
- Charge the unit at least every 3 month in event it is not used regularly.
- Please replace the pads after using 15-20 times (Varies widely depending on use, body type and care)
- Do not apply pads to the same area more than 30 minutes at a time.

Recommend massage points

Troubleshooting

Problem	Cause	Solution
During the massage, the skin feels a painful burning sensation or the stimulation becomes weakened.	Adhesive gel pads are not adhering firmly to the skin. Or the gel pads are too dry.	Apply a few drops of water to adhesive surface of each pad and make sure the pads are pressed firmly to the skin during application.
Unit is on, but no sensation is felt in the pads.	Both pads are not firmly on the skin. Or there is a loose connection. Or the intensity level needs to be increased.	Make sure both pads are firmly pressed to the skin. Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.
Adhesive gel pads do not stick to skin even after cleaning and moistening the gel pad.	Adhesive gel pads need to be replaced.	Replace the set of pads. Pads can be purchased from local distributor.
Unit does not turn on or the screen is dim.	Battery is low.	Recharge the battery.
Does not provide adequate pain relief.	Not using it long enough. Not lying down then massaging neck.	Use your massager 20-30 minutes at a time, 3-6 times per day. Lie on your back for neck support when you need relief for headaches or neck pain.

Manufactured by:

Shenzhen Hongqiangxing Electronics Co., Ltd.
Address: 4/F, Jingcheng Building, Xicheng
Industrial District, Xixiang Road, Baoan, Shenzhen,
China

Distributed by:

Wellkang Ltd
ADD: Suite B, 29 Harley Street, London, W1G9QR,
United Kingdom
Support@belifu.co

Made in China

This device complies with part 15 of the FCC rules

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