

Hand Massager

USER MANUAL



DARTWOOD

a WASSERSTEIN brand

About the Dartwood Hand Massage


This hand massager provides a deep-tissue massage for hands, wrists and fingers through multiple massage nodes and pads. It includes a heat function for added relaxation, improved circulation, reduced stress, and increased flexibility.

CAUTION

Always consult your doctor before using this hand massager or any vibration device.

May not be suitable if you recently had hand surgery, swelling, inflammation, congestion, severe arthritis, trauma, or other hand diseases. Do not use if you have had heart, brain, or lung disease, blood diseases, bleeding tendency, osteoporosis, soft tissue diseases or stroke. Never use on infants or if you are pregnant.

In the Box



1 x hand massager




1 x charging cable
(adaptor not included)

Safety Information

Read this manual carefully and keep it for future reference.

- For adult use only. Do not use if injured.
- Remove all the foreign objects like rings and bracelets from your hand before use, sharp objects will damage the airbag inside.
- Never leave the massager operating or charging unattended.
- Do not immerse in water. Keep away from liquids or heat sources.
- Do not use while charging.
- Each treatment should not be more than 30 minutes. Let the massager cool down for 10 minutes before reusing it.
- Do not drop.
- Do not use when driving.
- Do not attempt to repair, disassemble, or tamper with in any way.
- Only recharge with the original cable provided.
- Check the device regularly. Do not use if damaged.

Introduction to Parts




Before You Use

- Remove any jewelry or accessories from your hands (watch, ring, bracelet, etc.).
- Always start with lower intensity and find the most comfortable pressure for you.
- If you feel pain or discomfort, stop the massager immediately by pressing the power button.

How to Charge


1. Fully charge the device before use. Do not use it while charging.
2. Connect the charging cable to the USB-A socket (adapter not included) and the charging port at the side of the device.




3. Charging time is approximately 1.5 hours.
4. The power button light will pulsate when charging, and turn to solid when the battery is fully charged.
5. Remove the charging cable when fully charged.
6. Fully charge the battery before storing it for long periods of time to optimize the battery's healthiness.

How to Use

1. Place your hand inside the massager, fit your palm to the center of the first airbag. You can also adjust your preferred hand position.




2. Turn on the device by tapping the power button on the touch panel.




3. The display screen will turn on, indicating that it is standing by. The display will turn off if no function is selected within 5 seconds.
4. Select your hand size by tapping the hand icon on the device and start the treatment.

5. The massager will start and the intensity light will turn green. The massage cycle is approximately 15 minutes.
6. Select your preferred intensity by tapping the intensity icon. A green light indicates low intensity, a red light indicates medium intensity, and a blue light indicates high intensity.



7. Turn on the heating function by tapping the heating icon, the brightness indicates the level of heat. Do not use the heating for a long period of time to prevent low-temperature burn.



8. The massager will automatically turn off once the massage cycle has finished. You may also tap the power button again to turn the device off at anytime.

Frequently Asked Questions

1. **What are the benefits of this hand massager?**

This hand massager can provide a deep-tissue massage for hands, wrists and fingers through multiple massage nodes and pads. It can improve circulation, reduce stress and increase flexibility.

2. **My hand massager does not turn on, why?**

Make sure your massager is fully charged, connect the massager to a USB-A power source. The power button light will pulsate when charging, and turn to solid when the battery is fully charged.

3. **My fingers feel painful or get bruised when using the massager, what should I do?**

Make sure to place your hand correctly, your palm should sit on the first airbag and all your fingers should stay inside the massager around the second airbag. If you still find it is painful, please contact our customer service for further assistance.

Specifications

Model:	DW-HM908
Battery capacity:	1800 mAh
Massage method:	compression and heat
Cycle time:	15 min
Battery life:	7 - 9 cycles **
Charging time:	90 min
Charging input:	DC 5 V 2 A
Color:	white
Inner material:	antibacterial and anti-odor fabric
Product size:	8.3 x 7.0 x 4.1 in
Net weight:	1.6 lb

* Specifications are subject to change without notice.

** Battery life is subject to the intensity and actual environmental factors.

Get More Support

Scan the QR code below or email us at
contact@wasserstein-home.com



CONTACT US



DARTWOOD

a WASSERSTEIN brand

www.wasserstein-home.com

Made in China