Premium Automatic

Blood Pressure Monitor

INSTRUCTION MANUAL

Model BP3MS1-3B

Customer Help Line
1-800-568-4147

Automatic Blood Pressure Monitor - Instruction Manual

Table of contents

1. Introduction

- 1.1. Features of your blood pressure monitor
- 1.2. Gentle+™ Measurement System
- 1.3. Important information about self measurement

2. Important Information on Blood Pressure and its Measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done if regular high/low values are obtained?

3. Components of your Blood Pressure Monitor

4. Using your Monitor for the First Time

- 4.1. Installing the batteries
- 4.2. Setting the date and time
- 4.3. Cuff tube connection
- 4.4. Select the user

5. Measurement Procedure

- 5.1. Before measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Irregular Heartbeat Detector
- 5.6. Hypertension Indicator
- 5.7. Memory - Displaying the last 99 measurements
- 5.8. Discontinuing a measurement
- 5.9. Battery Charge Indicator
- 5.10. Using the AC Adapter

6. Bluetooth® functions & application (App)

- 6.1. Bluetooth® Icon Indicator
- 6.2. iOS devices
- 6.3. Android OS Phones

7. Error Messages/Troubleshooting

8. Care and Maintenance

9. Limited Warranty

10. Standards

11. Technical Specifications

12. Contacts for Support

1. Introduction

1.1. Features of your Automatic Blood Pressure Monitor

Your blood pressure monitor is an automatic digital blood pressure measuring device for use by adults on the upper arm at home or in your doctor's/nurse's office. It enables very fast and reliable measurement of systolic and diastolic blood pressure as well as pulse through the oscillometric method. This device offers clinically proven accuracy and has been designed to be user friendly.

Before using, please read this instruction manual carefully and then keep it in a safe place. Please contact your doctor for additional questions on the subject of blood pressure and its measurement.

1.2. Gentle+™ Measurement System

Gentle+™, our most advanced system to date, combines clinical accuracy and premium features including:

- A quieter, more comfortable measurement experience that reduces the amount of time your arm is under pressure. Less pressure = more comfort!
- An ergonomic cuff design that is easier to use and fits a wider range of arm sizes (for arms 8.7"-18.1").

Attention!

1.3. Important information about self-measurement

- **Self-measurement means control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers.
- In cases of irregular heartbeat, measurements made with this instrument should only be evaluated after consultation with your doctor.

Electromagnetic interference:

The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

2. Important Information on Blood Pressure and its Measurement

2.1. How does high or low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and speed of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is controlled by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the “blood ejection” (Systole) the value is highest (systolic blood pressure value). At the end of the heart’s “rest period” (Diastole) pressure is lowest (diastolic blood pressure value).

Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is very high if your diastolic pressure is above 100 mmHg and/or your systolic blood pressure is over 160 mmHg, **while at rest**. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

If your systolic blood pressure values are between 140 mmHg and 159mmHg and/or the diastolic blood pressure values between 90 mmHg and 99mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**

Which values are normal?

The following table is for classifying blood pressure values in adults in accordance with the National Institutes of Health JNC7 (Data in mmHg):+

Category	Systolic (mmHg)	Diastolic (mmHg)
Normal	<120	and <80
Pre-Hypertension	120-139	or 80-89
Hypertension		
Stage 1 Hypertension	140-159	or 90-99
Stage 2 Hypertension	≥160	or ≥100

Additional information

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension.” Consult your doctor.
- **Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.**
- **This device is not intended for use on children 12 years of age or younger.** Please ensure that children are supervised around this device, some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

2.3. What can be done if regular high or low values are obtained?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs

(heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.

- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

A) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Consistently follow all medical instructions for treating illness such as:
- Diabetes (Diabetes mellitus or sugar diabetes)
- Fat metabolism disorder
- Gout

C) Habits


- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea, chocolate).

D) Your physical condition

- **After a medical examination, and with your doctor's approval and direction; exercise.**
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. Speak with your doctor regarding the type and extent of exercise that is appropriate for you.

3. Components of your blood pressure monitor


a) Measuring unit

b) One Size Easy Fit Cuff:

For arm circumference 22cm - 46 cm (8.7" - 18.1").

If you ever need to buy a replacement cuff, call 1-800-568-4147 for information.

Please Note:


Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call 1-800-568-4147 for additional information. You may prefer to use one of our wrist type monitors.

4. Using your Monitor for the First Time

4.1 Installing the batteries

Insert batteries immediately after unpacking the unit. The battery compartment is located on the bottom of the device.

- a) Remove cover.
- b) Insert the batteries (4 x size AA 1,5V), observing the indicated polarity.
- c) If a battery warning appears in the display, the batteries are discharged and must be replaced.


Low Battery
Indicator

Attention!

- After the low battery indicator appears, the device won't function until the batteries have been replaced.
- Please use "AA" Long-Life or Alkaline 1.5 V batteries and replace them all at once.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.

4.2 Setting the date and time


1. After the new batteries are activated, the year number flashes in the display. You can advance the year by pressing the "M" button. To confirm and then set the month, press the time button.
2. You can now set the month using the "M" button. To confirm and then set the day, press the time button.
3. Please follow the instructions above to set the day, hour and minutes.
4. Once you have set the last minute and pressed the time button, the date and time are set and the time is displayed.
5. If you want to change the date and time, press and hold the time button down for approx. 3 seconds. Release and press the time button once more and the year will begin flashing. Now you can enter the new values as described above.

Memory
Button

Time/Dat
Button



4.3. Cuff tube connection

Insert the blue cuff tube into the blue opening on the side of the monitor indicated by the drawing of a cuff.

4.4. Select the User

This blood pressure monitor is designed to store 99 measurements for each of two users. Before taking a measurement, be certain that the correct user has been selected.

- a) With the unit off, press and hold the TIME button for 3 seconds until the user icon in the bottom right-hand corner of the LCD screen flashes.
- b) Press the "M" button to toggle between users.
- c) Press the power button to make your selection.

5. Measurement Procedure

Please note: You should always be seated before and during measurement.

5.1. Before measurement:

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about 5 minutes before taking a measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always compare measurements taken at the same time of day, since blood pressure changes during the course of the day, as much as 20-40 mmHg.

5.2. Common sources of error:

Note: Comparable blood pressure measurements always require the same conditions! Conditions should always be quiet.

- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured! Each 15cm (6 inch) difference in height between your heart and the cuff results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please call 1-800-568-4147.
- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results.

Consecutive blood pressure measurements should be repeated after a 1 minute pause or after your arm has been held up in order to allow the accumulated blood to flow away.


5.3. Fitting the Cuff

- a) The cuff is pre-formed for easier use. Remove tight or bulky clothing from your upper arm.
- b) Wrap the cuff around your upper left arm. The rubber tube should be on the inside of your arm extending downward to your hand. Make certain the cuff lies approximately 1/2" (1 to 2 cm) above the elbow.

Important! The red strip on the edge of the cuff (Artery Mark) must lie over the artery which runs down the inner side of the arm.

TIP: Align red artery mark to pinky finger.

- c) To secure the cuff, wrap it around your arm and press the velcro together.
- d) There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit.
- e) Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.
- f) Remain seated quietly for at least five minutes before you begin the measurement.


Comment:

If it is not possible to fit the cuff to your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Comparable blood pressure measurements always require the same conditions (Relax for 5 minutes before taking a measurement).

Please consult your physician for any of the following situations:

- The application of the cuff over a wound or inflamed area.
- The application of the cuff on any limb with intravascular access or therapy, or an arterio-venous (A-V) shunt.
- The application of the cuff on the arm on the side of a mastectomy.
- Simultaneous use with other medical monitoring equipment on the same limb.
- The blood circulation of the user needs to be checked.

Cuff on
right arm

5.4. Measuring procedure

After the cuff has been appropriately positioned on the arm and connected to the blood pressure monitor, the measurement can begin. Remain still during your measurement, do not flex muscles.

- a) Remain seated quietly for 5 minutes before you begin the measurement.
- b) Press the power button. The pump begins to inflate the cuff. In the display, the increasing cuff pressure is continually displayed.
- c) When the device has detected your pulse, the heart symbol in the display begins to blink.
- d) After reaching the inflation pressure, the pump stops and the pressure will release. The cuff pressure (large digits) is displayed during the measurement.
- e) When the measurement has been concluded, the measured systolic and diastolic blood pressure values, and the pulse are displayed.
- f) The appearance of the irregular heartbeat symbol signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements. **NOTE:** We recommend contacting your physician if you see this indicator frequently.
- g) The measurement results are displayed until you switch the device off. If no button is pressed for 1 minute, the device switches off automatically.

Pumping
Pressure

62

Measuring

46
♥

Systolic

Diastolic


Pulse

129
78
15

Measurement
complete

5.5. Irregular Heartbeat Detector

This symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure—repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily), we advise you to tell your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Irregular Heartbeat Symbol

This instrument is an oscillometric blood pressure monitor device that also analyzes pulse frequency during measurement. The instrument is clinically tested.


If pulse irregularities occur during measurement, the irregular heartbeat symbol is displayed after the measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

5.6. Hypertension Indicator


The colored bar on the left-hand edge of the display is the hypertension range. The black triangle on the left hand side of the screen will shift depending upon whether the readout value is within the normal (green), borderline (yellow) or danger (red) range.

The classification is based on standards established by the National Institutes of Health JNC7, 2003. Refer to the chart in section 2.2 of this instruction manual for details of the classifications.

With each reading, a corresponding black triangle is displayed to the right of the hypertension indicator.

- If your triangle is in the green zone, your measurement is "Normal," according to NIH standards.
- If your triangle is in the yellow zone, your measurement is "Pre-Hypertension."
- If your triangle is in the orange zone, it is "Stage 1 Hypertension."
- If your triangle is in the red zone, it is "Stage 2 Hypertension."

5.7. Memory – Displaying the last 99 Measurements

At the end of a measurement, this monitor automatically stores each result with date and time. This unit stores 99 memories for each of 2 users.

Viewing the stored values

With the unit off, press the "M" button. The display first shows "A", then shows an average of all measurements stored in the unit. Please note: Measurements for each user are averaged and stored separately. Be certain that you are viewing the measurements for the correct user.

Pressing the "M" button again displays the previous value. To view a particular stored memory, press and hold the "M" button to scroll to that stored reading.

Memory full

When the memory has stored 99 results for the selected user, a new measured value is stored by overwriting the oldest value.


Clear all values

If you are sure that you want to permanently remove all stored values, hold down the "M" button (the instrument must have been switched off beforehand) until "CL" appears and then release the button. If you do not want to clear the values, press the START/STOP button. To permanently clear the memory, press the "M" button while "CL" is flashing.

Individual values cannot be cleared.

5.8. Discontinuing a Measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the power button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

5.9. Battery Charge Indicator

Batteries almost discharged


When the batteries are approximately 75% used the battery symbol will flash as soon as the instrument is switched on (if at least one of the batteries still has some charge).

Although the instrument will continue to measure reliably, you should obtain replacement batteries.

Batteries discharged— replacements required

When the batteries are discharged, the battery symbol will appear unblinking as soon as the instrument is switched on. You cannot take any additional measurements and must replace the batteries.

1. Open the battery compartment at the back of the instrument by pushing inward at the two arrows and pulling it out.
2. Replace the batteries at the same time – ensure they are correctly connected, as shown on the symbols in the compartment.
3. The memory retains all values although date and time must be reset - the year number therefore flashes automatically after the batteries are replaced.
4. To set date and time, follow the procedure described in Section 4.2.

Which batteries and which procedure?

Use four new, longlife 1.5V AA batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period the batteries should be removed.

5.10. Using the AC Adapter


You may also operate this monitor using the included AC Adapter (output 6V DC/600 mA with DIN plug).

Use only the included AC Adapter to avoid damaging the unit.

1. Ensure that the AC adapter and cable are not damaged.
2. Plug the adapter cable into the AC Adapter port on the blood pressure monitor.
3. Plug the adapter into your electrical outlet.

When the AC adapter is connected, no battery current is consumed.

Note: No power is taken from the batteries while the AC Adapter is connected to the monitor. If electrical power is interrupted, (e.g., by accidental removal of the AC adapter from the outlet) the monitor must be reset by removing the plug from the socket and reinserting the AC adapter connection. Please consult us if you have any questions regarding the AC adapter by calling 1-800-568-4147.

6. Bluetooth® functions and application (App)

Please download the Microlife Connected Health App

(referred to as App in this document) from Apple's App Store® or Google Play® before pairing your devices. The App can be executed on a number of mobile platforms with the following specifications:

Compatibility:

- Apple® devices running iOS 8.0 or higher on iPhone®, iPad®, and iPod touch®.
- Android® phones running Android® 4.4.2 or higher.

Questions?

- To access the tutorial, go to the App home screen (Microlife Connected Health) and press the "Settings" button.
- Press "Tutorial" and select your topic.
- If you have any issues with your setup, please contact customer service toll free: 1-800-568-4147 for additional support

6.1 Bluetooth® Icon Indicator

The bluetooth® icon indicator on your blood pressure monitor, located in the upper left portion of the screen, is designed to provide information about the connection between your device and the monitor. Please note the following:

- The bluetooth® icon on your blood pressure monitor will flash initially when ready to pair with a device, and continue flashing for 2 minutes until pairing occurs,
- If pairing does not occur within 2 minutes, the Bluetooth® function on the monitor will automatically turn off.
- If the bluetooth® icon flashes for 10 seconds and then turns off, there is a problem with the connection.
- The icon will remain unblinking on the monitor when data is uploading to the App.

Bluetooth
Icon

Troubleshooting: Bluetooth® connection:

- Monitor: Turn off and wait for 1 minute.
- Device (eg. phone): Return to the Main Menu and reselect Download Data again.

If connection issues persist, please contact customer support toll free: 1-800-568-4147.

6.2 iOS devices:

Setting up your iOS device:

- In the settings menu, turn Bluetooth® “ON”
- Select the Microlife Connected Health App
- Select “User Profiles”
- Select “+” at the top of the screen to create a new user profile.
- Review the terms of service, and select “Accept” to proceed.
- Enter your information and select the check mark to confirm.


Pairing the devices for the first time:

On your device:

- Open the Microlife Connected Health App and select “Settings.”
- Select “Pair Device.”

On your blood pressure monitor:

- Press and hold the power button on your monitor for 5 seconds. Once connected, a Bluetooth® icon will appear on the blood pressure monitor’s screen.
- “OK” will be displayed if the pairing is successful.

Bluetooth icon

Setting up your blood pressure monitor:

Confirm the correct user is selected on your blood pressure monitor:


- Confirm the user icon on the bottom right hand side of the monitor screen is set to user 1 or 2.
- To change users on your monitor, see section 4.4.

Saving a measurement automatically (iOS):

- Open the Microlife Connected Health App
- Select "Download data."
- Take your blood pressure.
- Once complete, the blood pressure monitor will automatically activate a Bluetooth® connection with your iOS device. A Bluetooth® icon will appear on the monitor's screen when connected. (The devices must have been previously paired, and be in range).
- An image of a blood pressure monitor and a phone connecting will appear, and the readings will transfer.
- To save your data and return to the main menu, select the check mark in the upper right hand corner.

Your measurements will now be available for review through the "View Data" button on the home screen.

Bluetooth icon

Data Download

Saving measurements manually (iOS):

1. Open the Microlife Connected Health App
2. Select "User Profiles," and select the profile which you would like to upload. Once selected, the profile will be highlighted, and a check mark will appear to its left.
3. Use the home button to exit and return to the main menu.
4. Select "Download data." An image of a blood pressure monitor and a phone connecting will appear, and the readings will transfer.
5. To save your data and return to the main menu, select the check mark in the upper right hand corner.

Your measurements will now be available for review through the "View Data" button on the home screen.

User Profiles

Download
data

Saving readings stored on your monitor (iOS):

- Open the Microlife Connected Health App
- Select the profile which you would like to upload. Once selected, the profile will be highlighted, and check mark will appear to its left.
- On your blood pressure monitor: Press and hold the power button for 5 seconds. Once connected, a Bluetooth® icon will appear on the blood pressure monitor's screen.
- Follow steps 4 & 5 in "saving measurements manually"

6.3 Android OS Phones

Setting up your phone:

- In the settings menu, turn Bluetooth® "ON"
- Select the Microlife Connected Health App
- Select "User Profile"
- Select "+" at the top of the screen to create a new user profile.
- Review the terms of service, and select "Accept" to proceed.
- Enter your information and select the check mark to confirm.

On your blood pressure monitor:

- Confirm the correct user is selected on your blood pressure monitor by checking the icon on the bottom right hand side of the monitor screen is set to your preference (user 1 or 2).
- To change users on your monitor, see section 4.4 in this manual.

Saving readings to your Android phone:


On your phone:

- Open the Microlife Connected Health App
- Select "User Profiles," and then select the profile which you would like to upload (the profile will be highlighted in gray when selected).

On your blood pressure monitor:

- Press and hold the power button on your monitor for 5 seconds. Once connected, a Bluetooth® icon will appear on the blood pressure monitor's screen.

Bluetooth icon

On your phone:

- Select "Download data." Your phone will search and list available Bluetooth® devices.
- Select your monitor from the list (BP3MS1-3B), and the readings will automatically transfer.
- After the readings have been transferred, a summary of all your readings will be displayed.
- To save your readings and return to the main menu, select the check mark in the upper right hand corner.

Download data

Note: If no connection between the monitor and the phone has been made within 10 seconds, a "scan" button will appear in the top right corner of the phone's screen

- To search for your monitor again, press "scan."
- Once the phone has detected your monitor, select it from the list (BP3MS1-3B).

7. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.

(Example Error no. 1)

Error No.		Possible cause(s) / Solutions
ERR 1		The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. This situation can occur, for example, if the tube becomes unplugged after the systolic blood pressure has been measured. Additional possible cause: No pulse has been detected.
ERR 2		Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artifact).
ERR 3		Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.
ERR 5		The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions. Contact your doctor if this message occurs frequently.
HI		The cuff pressure is too high. Relax for 5 minutes and repeat the measurement.*
LO		The pulse is too low (less than 40). Repeat the measurement.*

***If this or any other problem occurs repeatedly, please consult your doctor.**





Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction	Remedy
The display remains blank when the instrument is switched ON although the batteries are in place.	• Check battery installation. • If the display is unusual, remove the batteries and then exchange them for new ones • Check polarity
The pressure does not rise although the pump is running.	• Check the connection of the cuff tube and connect properly.
The device frequently fails to measure the blood pressure values, or the values measured are too low or too high.	1. Fit the cuff correctly on the arm. (see 5.3) 2. Before starting the measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take off articles of clothing if necessary. 3. Measure blood pressure again in complete peace and quiet.
Every measurement results in different values, although the device functions normally and normal values displayed.	• Please read the following information and the points listed under "Common sources of error".(5.2) Repeat the measurement. Blood Pressure changes constantly. The observed readings may accurately reflect your pressure.
The blood pressure values differ from those measured by my doctor.	• Record daily values, consult your doctor. • Pressure readings in your doctor's office may be higher due to anxiety. • Your blood pressure fluctuates throughout the day so record the time along with your reading.
After the instrument has inflated the cuff the pressure falls very slowly or not at all. (No reasonable measurement possible).	• Check cuff connections. Ensure the unit has not been tampered with.

8. Care and Maintenance

- a) Do not store the unit under direct sunlight, at a high temperature, or in high humidity or dust. Performance may be degraded.
- b) The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **The cuff with bladder must not be washed in a dishwasher, clothes washer, or submerged in water.**
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the monitor!** This invalidates the manufacturer's warranty.
- g) Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Additional information

Blood pressure is subject to fluctuations even in healthy people.

Comparable measurements always require the same conditions (quiet conditions)! If fluctuations in readings are larger than 15 mmHg, consult your doctor.

Never attempt to repair the instrument yourself. Any unauthorized opening of the instrument invalidates all warranty claims.

Do not use the device if you think it is damaged or if anything appears unusual.

9. Limited Warranty

Your Blood Pressure Monitor is warranted for 5 years by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase.

The 5 year warranty applies to the monitor only. The following accessories are warranted for 1 year: cuff and adapter. Batteries are not covered by this warranty. There are no user serviceable parts inside.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents. Professional use, not following the operating instructions, and alterations made to the monitor or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A monitor or accessory as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A monitor or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost.

Please use the below customer service contact information to reach Microlife USA Inc. regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

10. Standards

Device standard:

Device corresponds to the requirements of the standard for non-invasive blood pressure monitors:

AAMI / ANSI / IEC 80601-2-30

IEC 60601-1, IEC 60601-1-2, IEC 60601-11

Bluetooth® system to V4.0, FCC Part 15C

Electromagnetic compatibility:

Device fulfills the stipulations of the International standard IEC 60601-1-2

IP20 Classification:

Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Keep Dry. 

Clinical testing:

Testing per the ESH Protocol and ISO 81060-2:2013 was conducted on blood pressure units using the same measurement technology.

Expected Service Life:

Monitor: 5 years

Cuff: 2 years.

Do NOT use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a distance from such devices when using this unit.

Please note: According to international standards, your monitor should be checked for accuracy every 2 years.

Trademark Usage:

Apple, the Apple logo, iPad, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Android and Google Play are both trademarks of Google Inc.

The Bluetooth® Smart word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks in this Blood Pressure Monitor is under license. Other trademarks and trade names are those of their respective owners.

©2016 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, LEADER, and the Leader LOGO are trademarks or registered trademarks of Cardinal Health. All other marks are property of their respective owners.

FCC Compliance:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications to the product are not approved by Microlife USA and could void the user's authority to operate the equipment under FCC jurisdiction.

This equipment has been verified to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by increasing the distance between the product and the affected device; or if applicable, reorientation the receiving antenna, or moving the product's power plug to another receptacle.

Models with wireless comply with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Federal Communications Commission (FCC) Statement

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) this device may not cause interference and
- 2) this device must accept any interference, including interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

11. Technical specifications:

Weight:	400 g (with batteries)
Size:	3.5" (W) x 6" (L) x 1.5" (H)
Storage temperature:	-20 to +55°C (-4°F to +131°F)
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C (50°F - 104°F)
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS:	60 to 255 mmHg
DIA:	40 to 200 mmHg
Pulse:	40 to 199 per minute
Cuff pressure display range:	0–299 mmHg
Memory:	Automatically stores the last 99 measurements for 2 users (total 198)
Measuring resolution:	1 mmHg
Accuracy:	pressure within ± 3 mmHg or 2% of reading > 200 mmHg pulse ± 5 % of the reading
Power source:	a) 4 AA batteries, 1.5 V b) AC adapter 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
Accessories:	One size easy fit cuff 22-46 cm (8.7" - 18.1") Technical alterations reserved!

Made in China

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.

Type BF applied part.

***Not for use on children 12 years of age or younger.**

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Android y Google Play son marcas comerciales de Google Inc.